


# Assay Method Using Peptidyl-pNA Substrates

## Principle

A protease with limited specificity hydrolyzes a peptidyl-pNA substrate, releasing *p*-nitroaniline (*p*NA) as follows:



The initial rate of increase in the *p*NA concentration can be monitored photometrically at 405 nm.




Fig. UV-Absorption Spectra


## Reagents

- 1) Substrate stock solution: 10 mM (DMSO or distilled water)
- 2) Buffer
- 3) Enzyme solution

## Procedure

Choose the proper conditions for the measurement, such as substrate concentration and sensitivity setting, depending on the purpose of the experiment and the instrument available. Described here is one of the recommended procedures for the photometric method (initial-rate method).

- 1) Set a spectrophotometer at  $\lambda = 405 \text{ nm}$  at  $25^\circ \text{C}$  ( $\epsilon_{405 \text{ nm}} = 9920$ )\*
- 2) Pipette 2940  $\mu\text{l}$  of buffer and 30  $\mu\text{l}$  of substrate stock solution into the cuvette
- 3) Incubate in the spectrophotometer for 3-4 min (for temperature equilibration)
- 4) Add 30  $\mu\text{l}$  of enzyme solution
- 5) Record the increase of the absorption intensity for 3-4 min
- 6) Calculate the amount of *p*NA released using the following equation



\*R. Lottenberg, U. Christensen, C.M. Jackson, and P.L. Coleman, *In*, Proteolytic Enzymes Part C, *Methods in Enzymology*, Vol. 80, (L. Lorand, ed.) Academic Press, New York, 1981, pp. 341-361.