

## Replacement Screw Caps and Septa

CODE	PRODUCT	USD
------	---------	-----

### For Manual Reaction Vessels

CAP-22201-PI	5,10, and 15 mL	4
CAP-22203-PI	30 mL - 400 mL	6
CAP-22205-PI	600 mL and Larger	8



### Replacement Screw Cap (Open) Plus Teflon-Coated Septum

CAP-22225-PI	10 and 15 mL	4
CAP-22232-PI	30 mL - 400 mL	6



### Extra Septa (Teflon-Coated Silicone Rubber)

SEP-22825-PI	10 and 15 mL	2
SEP-22832-PI	30 mL - 400 mL	2
SEP-22845-PI	600 mL and Larger	4

## Reaction Vessels for PS-3 Synthesizer

CODE	PRODUCT	USD
------	---------	-----

SHG-21003-PI		156
--------------	--	-----



## Evaporator Trap (Thick Walled) (Splash Guard Adapter)

SHG-20200-PI	200 mL	75
--------------	--------	----



## Buechner Funnel

GBF-20250-PI	70 mL Buechner Funnel Recommended for Filtration of CLEAR Resins (<2 grams)	77
--------------	---	----



GBF-20255-PI	125 mL Buechner Funnel Recommended for Filtration of CLEAR Resins (2-5 grams)	85
--------------	---	----

## Manual Solid-Phase Synthesis Shakers


CODE	PRODUCT	USD
------	---------	-----

SMS-22700-PI	120° Manual Synthesis Shaker	1029
SVR-22701-PI	180° (or 120°) Degree Variable Rate Shaker	1095
SVR-22702-PI	European Conversion Upgrade for Shaker	100



## Orbital Shaker with Universal Adapter

SHO-22500-PI	Orbital Shaker	895
SHO-22501-PI	Universal Adapter for Orbital Shaker	290




## Teflon Tubing Connectors

STC-22900-PI	(for attaching 1/4" glassware to 1/8" Teflon tubing)	49
STC-22902-PI	(for attaching 1/4" glassware to 1/4" Teflon tubing) **two connectors required per vessel	49



STB-21018-PI	Teflon Tubing 1/8" (10 feet)	25
--------------	------------------------------	----

The 1/8" diameter Teflon tubing is ideal for connecting a semi-automated synthesis assembly. Use this tubing with the available Teflon connectors to link 1/4" glass tubing (reaction vessels). Typically, solid phase assemblies connect to valves, stopcocks, aspirators, pumps, nitrogen tanks, and reagent bottles through a variety of configurations. Purchase this tubing in ample lengths and let your imagination and creativity be your guide.  
(Or call Peptides International for a few examples of potential applications).



## Hour-Glass Reaction Vessels with Teflon Plug

CODE	PRODUCT	USD
SHG-21028-PI	28 mL	159
SHG-21040-PI	40 mL	169
SHG-21060-PI	60 mL	179
SHG-21080-PI	80 mL	189
SHG-21100-PI	100 mL	199
SHG-21120-PI	120 mL	209
SHG-21250-PI	250 mL	219



## Hour-Glass Reaction Vessels with Screw Cap Closure

CODE	PRODUCT	USD
SHG-22028-PI	28 mL	129
SHG-22040-PI	40 mL	139
SHG-22060-PI	60 mL	149
SHG-22080-PI	80 mL	159
SHG-22100-PI	100 mL	169
SHG-22120-PI	120 mL	189
SHG-22250-PI	250 mL	229



## Manual Reaction Vessels

SHG-20205-PI	5 mL	110
SHG-20210-PI	10 mL	110
SHG-20215-PI	15 mL	110
SHG-20230-PI	30 mL	130
SHG-20260-PI	60 mL	150
SHG-20290-PI	90 mL	165



## Manual Reaction Vessels with Baffles for Better Mixing

SHG-20295-PI	200 mL	225
SHG-20300-PI	300 mL	245
SHG-20400-PI	400 mL	260
SHG-20600-PI	600 mL	295



## Manual Reaction Vessels with Jacket for Temp. Control

SHG-22215-PI	15 mL	195
SHG-22230-PI	30 mL	195
SHG-22260-PI	60 mL	235
SHG-22290-PI	90 mL	265



## Manual Reaction Vessels with 3-Way Teflon Stopcock with Two 1/4" Outlets

SHG-22315-PI	15 mL	145
SHG-22330-PI	30 mL	175
SHG-22360-PI	60 mL	195
SHG-22390-PI	90 mL	210



## Amber Manual Reaction Vessels

SHG-21115-PI	15 mL	155
SHG-21130-PI	30 mL	175
SHG-21160-PI	60 mL	205
SHG-21190-PI	90 mL	219



## Merrifield Reaction Vessels

SHG-20330-PI	30 mL	155
SHG-20340-PI	40 mL	165
SHG-20375-PI	75 mL	175



## Manual Reaction Vessels with Vacuum Adapter

SHG-22630-PI	30 mL Manual Reaction Vessel with 24/40 Vacuum Adapter Joint	155
SHG-22660-PI	60 mL Manual Reaction Vessel with 24/40 Vacuum Adapter Joint	190
SHG-22690-PI	90 mL Manual Reaction Vessel with 24/40 Vacuum Adapter Joint	205
SHG-22695-PI	200 mL Manual Reaction Vessel with 24/40 Vacuum Adapter Joint	240



## Ananth Reaction Vessels with Two 40 mL Chambers and Two Closures

SRV-25040-PI	Teflon closure	259
SRV-25240-PI	Screw cap closure	220

